

DO YOU TRUST YOUR PROCESSES ?

Put your trust in the SteriTrust Touchscreen Traceability System

Upgrade your Midmark steriliser with a reliable, wireless traceability system.

The touchscreen technology is easy-to-use and automatic cloud backup provides peace of mind.

Data is quickly transferred to patient files using the barcode scanner, reducing paperwork and increasing accuracy.

SteriTrust can be installed by your autoclave service technician and is compatible with all practice management software.

CONTACT US TODAY

sales@icona.net.au

www.icona.net.au

1300 442 662



SteriTrust

ICONA
DESIGNED WITH YOU IN MIND

SteriTrust Touchscreen Traceability System

- Quick and simple traceability
- Easy-to-use touchscreen technology
- Compatible with all practice management software
- Wireless connectivity
- Automatic cloud backup
- Reduces paper work
- Increases accuracy

AT END OF CYCLE

A user name is assigned to cycle records

A detailed barcode label is printed

Extra labels can be printed

Cycle data can easily be transferred to your patient's file using the barcode scanner

LABELS TYPICALLY SHOW

Practice name

Cycle number

Autoclave serial number

Date and time

Batch expiry date

Cycle type

Result

LOGGING

All cycle data is saved to on-board storage

Automatic cloud backup

Cycle data can be restored from cloud backup

Put your trust in the SteriTrust Touchscreen Traceability System



ICONA
DESIGNED WITH YOU IN MIND

SteriTrust

sales@icona.net.au

www.icona.net.au

1300 442 662